



## Continence Care Products



Making a difference in the journey of life.

# Table of Contents

<b>Secure Start<sup>SM</sup> Services</b> .....	2-3
Intermittent Catheters .....	4-16
<b>Infyna Chic<sup>TM</sup></b> Intermittent Catheters.....	4-5
<b>VaPro<sup>TM</sup></b> No Touch Intermittent Catheters.....	6-9
<b>Onli<sup>TM</sup></b> Intermittent Catheters .....	10-11
<b>Advance Plus<sup>TM</sup></b> Intermittent Catheter Systems.....	12-13
<b>Apogee Plus<sup>TM</sup></b> Intermittent Catheter Systems .....	14-16
Male External Catheters .....	17-18
<b>InView<sup>TM</sup></b> Silicone Male External Catheters .....	17
<b>Extended Wear<sup>TM</sup></b> Latex Male External Catheters.....	18
<b>Everyday<sup>TM</sup></b> Latex Male External Catheters.....	18
Leg Bag Systems .....	19-20
Vented Urinary Leg Bag Combination Packs.....	19
Urinary Leg Bag Combination Packs .....	19
Urinary Leg Bags .....	20
Leg Bag Accessories .....	21
Fabric Leg Bag Straps.....	21
Extension Tubings .....	21
Collection Systems.....	22-24
Bedside Drainage Collection System with Anti-Reflux Valve.....	22
Male Urinary Pouch External Collection Devices .....	23
Female Urinary Pouch External Collection Device .....	24
Latex Statement.....	25
HCPCS Codes Disclaimer .....	25
Index .....	28

# Hollister Secure Start<sup>SM</sup> Services

Personalized Support



## We're here to help.

If you have questions about finding the right intermittent catheter, sorting out insurance coverage, or connecting with community resources, we're here to help. You don't have to figure it out on your own.



**Healthcare Professionals, enroll your patient today!**

**[www.Hollister.com/CCEnroll](http://www.Hollister.com/CCEnroll) | 1.888.808.7456**



Hollister Secure Start services offer free support for catheter users for as long as you need it, regardless of the brand of products you use.

**Call Us Today! 1.888.808.7456**  
**[www.securestartservices.com](http://www.securestartservices.com)**

Hollister Secure Start<sup>SM</sup> services offer free customized support for those who are prescribed intermittent catheters for bladder management, regardless of the brand of product used, including:



**Finding the right products.**

Whether you are new to catheterization or an experienced user looking for options, we can help.



**Navigating insurance coverage.**

Based on the insurance plan, we can explain options and help determine whether any catheter product, regardless of brand, is covered as well as allowable quantities.



**Accessing educational and lifestyle information.**

We have long-standing relationships with a variety of organizations and can connect you with people, knowledge and tools to help support the life you want to live.



**Identifying product supplier options.**

Upon request, we can provide you with supplier options who accept your insurance and who can provide your preferred product.



## Thoughtfully refined with women in mind

Subtle enhancements, improved outcomes.

We're committed to improving the self-catheterization experience for the women who choose us, which is why we're introducing some subtle yet significant changes to our Infyna Chic™ intermittent catheters. Delivering even more comfort and confidence – for your peace of mind.

**Call to request a sample today!**

**1.888.808.7456**



Infyna Chic™

Intermittent Catheter



Making a difference in the journey of life.



Scan to  
Request a  
Free Sample

## Infyna Chic™ Intermittent Catheters

**81xx-30**

**Box of 30**

- ▶ Discreet to Carry with Unique Combination of Color, Design and Small Size
- ▶ Flip-top Cap on Catheter Case is Easy to Open with Just One Hand
- ▶ Reseal Case After Use and Dispose of Without Leaks or Odor
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter with Enhanced Lubrication
- ▶ Just Right Catheter Stiffness Helps Support Easy, Touch-Free Insertion
- ▶ Catheter Length Helps Give Confidence Bladder is Fully Drained

### HCPCS: A4351

Stock No	Length (inches)	Size	Description
8108-30	5.5	8 Fr	Straight
8110-30	5.5	10 Fr	Straight
8112-30	5.5	12 Fr	Straight
8114-30	5.5	14 Fr	Straight



With 100% No Touch Protection,  
whatever you touch won't touch  
VaPro™ catheters.

The protective tip and sleeve of VaPro™ catheters have been tested and have been shown to prevent the transmission of bacteria that may cause UTIs. Don't settle for anything less than 100% No Touch Protection.

**Call to request a sample today!**  
**1.888.808.7456**



**VaPro™**  
No Touch Intermittent Catheter



Making a difference in the journey of life.



Coudé Tip



Scan to  
Request a  
Free Sample

## VaPro Pocket™ Coudé No Touch Intermittent Catheters

**77xxx-30**

**Box of 30**

- ▶ Full Length No Touch Catheter in a Pocket-Size Package
- ▶ Protective Tip and Sleeve Allow for 100% No Touch Catheterization
- ▶ Coudé Tip
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter
- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex

**HCPCS: A4352**

Stock No	Length (inches)	Size	Description
77124-30	16	12Fr	Coudé
77144-30	16	14Fr	Coudé
77164-30	16	16Fr	Coudé



Coudé Tip



Scan to  
Request a  
Free Sample

## VaPro™ Coudé No Touch Intermittent Catheters

**73xxx**

**Box of 30**

- ▶ Protective Tip and Sleeve Allow for 100% No Touch Catheterization
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter
- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex
- ▶ Coudé Tip

**HCPCS: A4352**

Stock No	Length (inches)	Size	Description
73124	16	12 Fr	Coudé
73144	16	14 Fr	Coudé
73164	16	16 Fr	Coudé





Scan to  
Request a  
Free Sample

## VaPro Pocket™ No Touch Intermittent Catheters

70xxx-30

Box of 30

- ▶ Full length No Touch Catheter in a Pocket-Size Package
- ▶ Protective Tip and Sleeve Allow for 100% No Touch Catheterization
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter
- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4353

Stock No	Length (inches)	Size	Description
70104-30	16	10 Fr	Straight
70124-30	16	12 Fr	Straight
70144-30	16	14 Fr	Straight
70164-30	16	16 Fr	Straight



Scan to  
Request a  
Free Sample

## VaPro Plus Pocket™ No Touch Intermittent Catheters

71xxx-30

Box of 30

- ▶ Full Length No Touch Catheter in a Pocket-Size Package
- ▶ Protective Tip and Sleeve Allow for 100% No Touch Catheterization
- ▶ Integrated 1000 mL Collection Bag
- ▶ Smooth Eyelets
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4353

Stock No	Length (inches)	Size	Description
71082-30	8	8Fr	Straight
71102-30	8	10Fr	Straight
71122-30	8	12Fr	Straight
71142-30	8	14Fr	Straight
71084-30	16	8Fr	Straight
71104-30	16	10Fr	Straight
71124-30	16	12Fr	Straight
71144-30	16	14Fr	Straight
71164-30	16	16Fr	Straight



## VaPro Plus™ No Touch Intermittent Catheters

**74xxx**

**Box of 30**

- ▶ Protective Tip and Sleeve Allows for 100% No Touch Catheterization
- ▶ Smooth Eyelets
- ▶ Integrated 1000 mL Collection Bag
- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter
- ▶ Not Made with Natural Rubber Latex



Scan to  
Request a  
Free Sample

**HCPCS: A4353**

Stock No	Length (inches)	Size	Description
74122	8	12 Fr	Straight
74142	8	14 Fr	Straight
74124	16	12 Fr	Straight
74144	16	14 Fr	Straight




### DID YOU KNOW...

VaPro™ No Touch Intermittent Catheters offer the right balance of ease of use and protection:

- Prehydrated, hydrophilic and ready to use
- Protective sleeve allows for the catheter to be gripped anywhere and is proven to prevent the transmission of pathogens that may cause UTIs
- Protective tip helps shield the sterile catheter during insertion from bacteria located within the first 15mm of the distal urethra
- Packaging is designed with finger holes to be easy to open
- VaPro Plus™ and VaPro Plus Pocket™ catheters have an integrated 1000ml collection bag for catheterization at home or on the go
- Offered in pocket-size packaging for easy transportation and discreet out of home usage

# Onli™ Intermittent Catheters



Scan to  
Request a  
Free Sample

## Onli™ Intermittent Catheters

82xxx-30

Box of 30

- ▶ Ready-to-Use, Prehydrated Hydrophilic Catheter with Enhanced Lubrication
- ▶ Easy-to-Open, Peel-Apart Packaging with Two Finger Holes
- ▶ Ergonomic Sliding Gripper\* and Length of Flexible Funnel Help Provide Control
- ▶ Catheter has a Balance of Flexibility and Rigidity that Helps Promote Easy, Touch-Free Insertion
- ▶ Clear Tube Provides Increased Visibility
- ▶ Doesn't Stain Common Fabrics\*\* if Contents are Spilled

### HCPCS: A4351

Stock No	Length (inches)	Size	Description
82081-30	7	8 Fr	Straight
82101-30	7	10 Fr	Straight
82121-30	7	12 Fr	Straight
82141-30	7	14 Fr	Straight
82084-30	16	8 Fr	Straight
82104-30	16	10 Fr	Straight
82124-30	16	12 Fr	Straight
82144-30	16	14 Fr	Straight
82164-30	16	16 Fr	Straight

\*Sliding gripper is only on 16" catheter

\*\*Common fabrics include cotton, denim, linen and polyester



**Easier to open:**  
two finger holes and  
peel-apart packaging



**Enhanced lubrication:**  
more lubricious hydrophilic  
coating for easier insertion  
and withdrawal



**Clear catheter tube:**  
increased visibility and a  
more appealing look



## Inspired updates for easier use

**Small improvements that add up to big results.**

Introducing the improved Onli™ intermittent catheter. We've built on the features users love with some small but smart improvements that will make a big difference to their self-catheterization experience.

**Call to request a sample today!**  
**1.888.808.7456**



**Onli**™  
Intermittent Catheter



Making a difference in the journey of life.

# Advance Plus™ Intermittent Catheters



Scan to  
Request a  
Free Sample

## Advance Plus™ Intermittent Catheter Systems

94xxx/95xxx

Box of 100

- ▶ Gel Reservoir
- ▶ Protective Tip
- ▶ 1500 mL Collection Bag with Large Handle
- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4353

Stock No	Length (inches)	Size	Description
94064	16	6 Fr	Straight
94084	16	8 Fr	Straight
94104	16	10 Fr	Straight
94124	16	12 Fr	Straight
94144	16	14 Fr	Straight
94164	16	16 Fr	Straight
94184	16	18 Fr	Straight
95124	16	12 Fr	Coudé
95144	16	14 Fr	Coudé
95164	16	16 Fr	Coudé



Scan to  
Request a  
Free Sample

## Advance Plus Pocket™ Intermittent Catheter Systems

93xxx

Box of 100

- ▶ Full Length Advance Plus™ Intermittent Catheter System in a Pocket-Size Package

HCPCS: A4353

Stock No	Length (inches)	Size	Description
93104	16	10 Fr	Pocket
93124	16	12 Fr	Pocket
93144	16	14 Fr	Pocket
93164	16	16 Fr	Pocket



## Advance Plus™ Intermittent Catheter System Kits

96xxx/97xxx

Box of 100

- ▶ Advance Plus™ Catheter
- ▶ Kit Contains Gloves, Underpad Drape, and Benzalkonium Chloride Swabsticks
- ▶ Not Made with Natural Rubber Latex



Scan to  
Request a  
Free Sample

### HCPCS: A4353

Stock No	Length (inches)	Size	Description
96064	16	6 Fr	Straight
96084	16	8 Fr	Straight
96104	16	10 Fr	Straight
96124	16	12 Fr	Straight
96144	16	14 Fr	Straight
96164	16	16 Fr	Straight
96184	16	18 Fr	Straight
97124	16	12 Fr	Coudé
97144	16	14 Fr	Coudé
97164	16	16 Fr	Coudé



## Apogee IC™ Intermittent Catheters

**10xx/10xxx/11xxx**

**Box of 30**

- ▶ Smooth Eyelets
- ▶ Not Made with Natural Rubber Latex



Scan to  
Request a  
Free Sample

**HCPCS: A4351 for all except**

**HCPCS: A4352 for xxx26**

Stock No	Length (inches)	Size	Description
10806	6	8 Fr	Firm
11006	6	10 Fr	Firm
11206	6	12 Fr	Firm
11406	6	14 Fr	Firm
11446*	6	14 Fr	Firm
10610	10	6 Fr	Firm
10810	10	8 Fr	Firm
11010	10	10 Fr	Firm
10616	16	6 Fr	Firm
10816	16	8 Fr	Firm
11016	16	10 Fr	Firm
11216	16	12 Fr	Firm
11416	16	14 Fr	Firm
11616	16	16 Fr	Firm
11816	16	18 Fr	Firm
10826	16	8 Fr	Coudé
11026	16	10 Fr	Coudé
11226	16	12 Fr	Coudé
11426	16	14 Fr	Coudé
11626	16	16 Fr	Coudé
11826	16	18 Fr	Coudé
1061	16	14 Fr	Soft
1065	16	14 Fr	Firm, Curved Package

\*No funnel Qty: 30/box except curved catheters: 50/box



Scan to  
Request a  
Free Sample

## Apogee Plus™ Intermittent Catheter Systems

**Bxx/Bxxx**

**Box of 100**

- ▶ Protective Tip
- ▶ Pre-lubricated
- ▶ Smooth Eyelets
- ▶ 1500 mL Collection Bag
- ▶ Not Made with Natural Rubber Latex

### HCPCS: A4353

Stock No	Length (inches)	Size	Description
B8F	16	8 Fr	Firm
B10F	16	10 Fr	Firm
B12F	16	12 Fr	Firm
B14F	16	14 Fr	Firm
B16F	16	16 Fr	Firm
B18F	16	18 Fr	Firm
B12C	16	12 Fr	Coudé
B14C	16	14 Fr	Coudé
B16C	16	16 Fr	Coudé
B12S	16	12 Fr	Soft
B14S	16	14 Fr	Soft
B14F-Female	8	14 Fr	Firm



# Apogee™ Plus Intermittent Catheters



## Apogee Plus™ Intermittent Catheter System Kits

**BxxB/BxxxB**

**Box of 100**

- ▶ Apogee IC™ Intermittent Catheter
- ▶ Kit Contains Gloves, Underpad Drape, Gauze, and PVP Swabsticks
- ▶ Not Made with Natural Rubber Latex



Scan to  
Request a  
Free Sample

### HCPCS: A4353

Stock No	Length (inches)	Size	Description
B6FB	16	6 Fr	Firm
B8FB	16	8 Fr	Firm
B10FB	16	10 Fr	Firm
B12FB	16	12 Fr	Firm
B14FB	16	14 Fr	Firm
B16FB	16	16 Fr	Firm
B18FB	16	18 Fr	Firm
B12CB	16	12 Fr	Coudé
B14CB	16	14 Fr	Coudé
B16CB	16	16 Fr	Coudé
B12SB	16	12 Fr	Soft
B14SB	16	14 Fr	Soft
B14FB-Female	8	14 Fr	Firm

# Male External Catheters



Scan to  
Request a  
Free Sample



## InView™ Silicone Male External Catheters

**974xx/975xx/976xx**  
**974xx-100/975xx-100/976xx-100**

**Box of 30**  
**Box of 100**

- ▶ Self-Adhesive
- ▶ Clear Material Enables Monitoring of Skin
- ▶ Not Made with Natural Rubber Latex

**HCPCS: A4349**

**SPECIAL: shorter length sheath**

<b>Stock No Box of 30</b>	<b>Stock No Box of 100</b>	<b>Color</b>	<b>Size</b>
97425	97425-100	■	25 mm
97429	97429-100	■	29 mm
97432	97432-100	■	32 mm
97436	97436-100	■	36 mm
97441	97441-100	■	41 mm

**STANDARD: the daily catheter for general use**

<b>Stock No Box of 30</b>	<b>Stock No Box of 100</b>	<b>Color</b>	<b>Size</b>
97525	97525-100	■	25 mm
97529	97529-100	■	29 mm
97532	97532-100	■	32 mm
97536	97536-100	■	36 mm
97541	97541-100	■	41 mm

**EXTRA: 70% more adhesive than Standard for added security**

<b>Stock No Box of 30</b>	<b>Stock No Box of 100</b>	<b>Color</b>	<b>Size</b>
97625	97625-100	■	25 mm
97629	97629-100	■	29 mm
97632	97632-100	■	32 mm
97636	97636-100	■	36 mm
97641	97641-100	■	41 mm

# Male External Catheters



Scan to  
Request a  
Free Sample

## Extended Wear™ Latex Male External Catheters

920x

Box of 30

- ▶ Self-Adhesive
- ▶ Unique Inner Flap
- ▶ Easy Applicator Collar
- ▶ Contains Natural Rubber Latex

HCPCS: A4349

Stock No  
Box of 30

Color

Size

9206	■	22 mm-25 mm
9207	■	26 mm-30 mm
9208	■	31 mm-35 mm
9209	■	36 mm-39 mm

**Caution:** This product contains natural rubber latex which may cause allergic reactions.

## Everyday™ Latex Male External Catheters

910x

Box of 30

- ▶ Self-Adhesive
- ▶ Easy Applicator Collar
- ▶ Contains Natural Rubber Latex

HCPCS: A4349

Stock No  
Box of 30

Color

Size

9106	■	22 mm-25 mm
9107	■	26 mm-30 mm
9108	■	31 mm-35 mm
9109	■	36 mm-39 mm

**Caution:** This product contains natural rubber latex which may cause allergic reactions.



Scan to  
Request a  
Free Sample

# Leg Bag Systems



Scan to  
Request a  
Free Sample

## Vented Urinary Leg Bag Combination Packs

9645/9655

Box of 10

- ▶ Includes:
  - (1) Sterile, Vented Urinary Leg Bag — Medium or Large
  - (1) Sterile, Extension Tubing with Vacuum Relief Valve, (18" Long, Oval, 11/32" Effective Inside Diameter), with Connector
  - (1) Pair, Fabric Leg Bag Straps — Medium (23" Long, 1-1/8" Wide)
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4358

Stock No	Description	Length (inches)	Width
9645	Medium — 18 oz (540 mL)	10"	4-1/4"
9655	Large — 30 oz (900 mL)	11-1/8"	5"



Scan to  
Request a  
Free Sample

## Urinary Leg Bag Combination Packs

9348/9349

Box of 12

- ▶ Includes:
  - (1) Sterile, Urinary Leg Bag — Medium or Large
  - (1) Non-Sterile, Extension Tubing (18" Long, Oval, 11/32" Effective Inside Diameter) with Connector
  - (1) Pair, Fabric Leg Bag Straps — Medium (23" Long, 1-1/8" Wide)
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4358

Stock No	Description	Length (inches)	Width
9348	Medium — 18 oz (540 mL)	10"	4-1/4"
9349	Large — 30 oz (900 mL)	11-1/8"	5"



Scan to  
Request a  
Free Sample

## Urinary Leg Bags

**9814/9805**

**Box of 10**

- ▶ Includes:  
(1) Sterile, Urinary Leg Bag — Medium or Large
- ▶ Not Made with Natural Rubber Latex

**HCPCS: A9270**

**HCPCS: A4358 (if straps are added)**

Stock No	Description	Length (inches)	Width
9814	Medium — 18 oz (540 mL)	10"	4-1/4"
9805	Large — 30 oz (900 mL)	11-1/8"	5"



### DID YOU KNOW...

#### Hollister Leg Bags Offer Many Benefits:

- Comfort with soft, vinyl material
- Discreet profile with pleats
- Strap security with wide, woven fabric
- Minimal dexterity needed with easy-to-operate drain valve
- Reduced risk of pressure with protected strap connection buttons
- Helps decrease urine back flow with anti-reflux valve

# Leg Bag Accessories



Scan to  
Request a  
Free Sample

## Fabric Leg Bag Straps

9342/9343

10 Pair

- ▶ Not Made with Natural Rubber Latex

HCPCS: A5114

Stock No	Description	Length (inches)	Width
9342	Small	15"	1-1/8"
9343	Medium	23"	1-1/8"



Scan to  
Request a  
Free Sample

## Extension Tubings

9345/9346

Box of 10

- ▶ Kink-resistant with Connector  
Oval, 11/32" Effective Inside Diameter
- ▶ Thermoplastic Elastomer
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4331

Stock No	Description	Length (inches)
9345	Non Sterile	18"
9346	Sterile	18"



Scan to  
Request a  
Free Sample

## Bedside Drainage Collection System with Anti-Reflux Valve

**9839**

**Box of 20**

- ▶ Sample Port
- ▶ Sterile Fluid Pathway
- ▶ Not Made with Natural Rubber Latex

**HCPCS: A4357**

Stock No	Description	Length (inches)	Width
9839	68 oz (2000 mL)	12"	9-3/4"



## Male Urinary Pouch External Collection Devices\*

**9873**

**Box of 10**

- ▶ Cut-to-Fit, FlexTend™ Skin Barrier
- ▶ Thin and Flexible
- ▶ No Additional Tape Needed
- ▶ Odor-Barrier Film
- ▶ Convenience Drain
- ▶ Not Made with Natural Rubber Latex

### HCPCS: A4326

Stock No	Description	Length (inches)	Width
9873	Male Urinary Pouch	7-1/2"	4-1/2"



Scan to  
Request a  
Free Sample

**9811**

**Box of 10**

- ▶ Cut-to-Fit, SoftFlex™ Skin Barrier
- ▶ Odor-Barrier Film
- ▶ Convenience Drain
- ▶ Not Made with Natural Rubber Latex

### HCPCS: A4326

Stock No	Description	Length (inches)	Width
9811	Male Urinary Pouch	7-1/2"	4-1/2"

*\*Meant to be connected to leg bag or gravity drainage*





Scan to  
Request a  
Free Sample

## Female Urinary Pouch External Collection Device\*

**9840**

**Box of 10**

- ▶ Cut-to-Fit, SoftFlex™ Skin Barrier
- ▶ Odor-Barrier Film
- ▶ Convenience Drain
- ▶ Alternative to Indwelling Catheter for Bedridden Females
- ▶ Not Made with Natural Rubber Latex

### HCPCS: A4328

<b>Stock No</b>	<b>Description</b>	<b>Length (inches)</b>	<b>Width</b>
9840	Female Urinary Pouch	7-1/2"	4-1/2"

*\*Meant to be connected to leg bag or gravity drainage*

## HCPCS Codes Disclaimer

The reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Hollister does not guarantee coverage or payment for any products. The ultimate responsibility for proper coding, satisfying reimbursement requirements, and obtaining reimbursement remains with the provider. Coding and coverage policies and guidelines are complex, can vary from one carrier or region to another, and are updated frequently. Providers should check with their local carriers or intermediaries often and should consult with counsel, a reimbursement specialist, and/or the DMECS website for any coding, coverage, reimbursement, or billing questions.

*HCPCS codes shown in this Product Catalog are current as of March 2023.*

## Latex Statement

The specifications for the Hollister Continence Care products listed below include natural rubber latex as a component; our suppliers add natural rubber latex as part of the production process. If a Hollister product(s) found in this catalog is not listed below, the specification for that product(s) and packaging does not include natural rubber latex as a component, and our component suppliers do not add natural rubber as part of their production process.

## Current Stock Numbers:

9106	9206
9107	9207
9108	9208
9109	9209

Prior to use of ALL continence care products, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

All intermittent catheters are Rx Only.





# Index

<b>Stock No</b>	<b>Page</b>	<b>Stock No</b>	<b>Page</b>
10xx	14	95xxx	12
10xxx	14	96xxx	13
11xxx	14	9645	19
70xxx-30	8	9655	19
71xxx-30	8	97xxx	13
73xxx	7	974xx	17
74xxx-30	9	974xx-100	17
77xxx-30	7	975xx	17
81xx-30	5	975xx-100	17
82xxx-30	10	976xx	17
910x	18	976xx-100	17
920x	18	9805	20
93xxx	12	9811	23
9342	21	9814	20
9343	21	9839	22
9345	21	9840	24
9346	21	9873	23
9348	19	Bxx/Bxxx	15
9349	19	BxxB/BxxxB	16
94xxx	12		

## Contact Us!

Questions about our products or services?

**1.888.808.7456**

For orders only:

**1.800.323.4060**

**Fax 1.847.680.1017**

### **Hollister Incorporated**

2000 Hollister Drive  
Libertyville, Illinois 60048

**[www.hollister.com](http://www.hollister.com)**

Hollister, the Hollister logo, Advance Plus, Advance Plus Pocket, Apogee, Apogee IC, Apogee Plus, Everyday, Extended Wear, Flexend, Infyna Chic, InView, Onli, Secure Start, Softflex, VaPro, VaPro Plus, VaPro Plus Pocket and VaPro Pocket are trademarks of Hollister Incorporated.

© 2023 Hollister Incorporated. 922530-0423 US-01924

