


A Hollister Technology



Moderma Flex  
Soft Convex.  
A Flexible Fit.  
**Made Secure.**



Hollister Ostomy. **Details Matter.**




# Introducing the NEW Moderma Flex Soft Convex Range

## Designed to stay secure, no matter what shape the body takes

When it comes to finding the right solution for patients, fit comes first. Comfort and ease-of-use are also essential to help ensure positive outcomes - in the hospital and after patients leave your care. The Moderma Flex soft convex range has been expanded to deliver both convex support and flexibility to a wider group.

- The convex barrier helps create an optimal fit by providing gentle pressure around flushed or recessed stomas, while conforming to uneven body contours.
- Better conformity helps to ensure a closer seal and protect against leaks. Fewer leaks helps promote healthy peristomal skin.
- The viewing option allows you to check pouch position and observe output with the lift of a flap.
- The skin barrier has an adhesive border for enhanced security.
- Built-in belt tabs can help keep the pouch in place, for added security and peace of mind.
- Integrated soft convexity may minimise the need for additional accessories, and help simplify patient education and use.



**Moderma Flex soft convex** helps create an optimal fit by providing gentle pressure around the stoma, while conforming to body contours.

# Moderma Flex Soft Convex Drainable Pouch

**AF300 filter** 

Generous air flow, coupled with an odour-eliminating carbon element, help reduce pouch ballooning and odour.

**FlexWear skin barrier formulation** 

The FlexWear skin barrier is a good balance of strength and flexibility.

**Adhesive Border**

The skin barrier has an adhesive border for enhanced security.

**Viewing Option**

Allows you to check pouch position and observe output with the lift of a flap.

**Built-in belt tabs**

Can help keep the pouch in place, for added security and peace of mind.

**Soft and flexible convexity** 

Better conformity helps to ensure a closer seal and protect against leaks. Fewer leaks helps promote healthy peristomal skin.


**Lock 'n Roll Microseal Closure** 

The microseal closure contains over a thousand tiny microseals which securely engage when pressed together, creating a water, weight, and pressure-resistant seal.

**Choice of pre-cut and cut-to-fit skin barriers**

Offers the largest (55 mm) cutting area of any Hollister convex skin barrier. Designed to suit a wider range of stoma sizes.

# The Moderma Flex Soft Convex **Closed Pouch**





## Soft convexity may benefit people with stomas who:

- Have skin folds or scars around the stoma
- Have a retracted or flat stoma
- Are experiencing leakage with other pouches
- Have postoperative hernias



## Ordering Information

		Colour	Stoma Size	Pack Size	Order No		
<b>Moderma Flex Soft Convex Drainable Pouch with FlexWear Skin Barrier and Adhesive Border</b>							
	Midi Cut-to-Fit	Beige with Viewing Option	15-25mm	10	5681025		
			15-38mm	10	5681038		
			15-55mm	10	5681055		
	Midi Pre-Cut	Beige with Viewing Option	20mm	10	5681120		
			25mm	10	5681125		
			30mm	10	5681130		
			35mm	10	5681135		
			40mm	10	5681140		
			Maxi Cut-to-Fit	Transparent	15-25mm	10	5691025
					15-38mm	10	5691038
					15-55mm	10	5691055
			Maxi Pre-Cut	Transparent	20mm	10	5691120
25mm	10	5691125					
30mm	10	5691130					
35mm	10	5691135					
			40mm	10	5691140		
<b>Moderma Flex Soft Convex Closed Pouch with SoftFlex Skin Barrier and Viewing Option</b>							
	Midi Cut-to-Fit	Beige	15-25mm	30	2261025		
			15-38mm	30	2261038		
			15-55mm	30	2261055		
	Midi Pre-Cut	Beige	20mm	30	2261120		
			25mm	30	2261125		
			30mm	30	2261130		
			35mm	30	2261135		
			40mm	30	2261140		
	Maxi Cut-to-Fit	Beige	15-25 mm	30	2281025		
			15-38 mm	30	2281038		
			15-55 mm	30	2281055		
	Maxi Pre-Cut	Beige	20 mm	30	2281120		
			25 mm	30	2281125		
			30 mm	30	2281130		
			35 mm	30	2281135		
			40 mm	30	2281140		

To request FREE samples, please contact Hollister on **0800 521 377** or email [samples.uk@hollister.com](mailto:samples.uk@hollister.com) for further information.

**United Kingdom**  
Hollister Limited  
Rectory Court, 42 Broad Street  
Wokingham, Berkshire RG40 1AB

[www.hollister.co.uk](http://www.hollister.co.uk)